


Check In – Dosage Interval


This enables you to select a time between dosages to avoid overdosing. For example, paracetamol cannot be taken within 4 hours of previous the dose.



Click on the **Check In** icon.



Select the appropriate **resident**.



Select the medication you want to create a dosage interval for (e.g.: Paracetamol).




Dosage Interval

Click on the **Dosage Interval** icon.



Dosage Interval



A green tick and a box for you to insert the interval time will appear. The default time is set to 4 hours but if you wish to change this click in the box and do so. Continue check in as normal.

Should someone try to administer the drug within the interval time a warning message will appear on the screen.