

## Meds Round - Body Map

For administering Creams and other topical medications.




Click on the Meds Round icon.



Select the **Appropriate resident.** 

Cavilon durable barrier 3391G (28g) cream 2000

Select the medication to be administered.



The medication administration window will appear with the body map.

You can add a comment at this stage if you like.

www.emarx.co.uk 01442 850050





Click on the save icon.




## **Body Map**

Please press the + button, then click and drag the blue cross to the area of the body map where the medication has been applied.

Multiple areas can be added, if required.

ок

The Body Map instruction box appears. Click Ok.



Click the Plus icon and drag the cross to the appropriate location.



## Saving information

Body Map administration successfully saved.

ОК

Click save

www.emarx.co.uk 01442 850050





You can view the body map point history by clicking on this icon in the preparation window.

<u>www.emarx.co.uk</u> 01442 850050