


Meds Round – Scanning Issues


In the case of a scanner malfunction or a damaged label, barcodes can be manually entered to confirm the correct medication. All broken scanners should be reported to the eMARx support team immediately. If the label is completely unreadable then the Override Scanning option can be used.



Click on the Meds Round icon.



Located next to the barcode will be a list of characters. This is what the barcode represents.



The barcode text can be typed directly into the Scan Medication Barcode window. Once entered click "Search".



If the barcode and number are unrecognisable then mark tick the override scanning icon in the check-in screen.