


## Stock Control – Toggle Automatic Stock Carry Forward On/Off

We'd recommend this method of carrying stock levels from one cycle to another. However, there are two other methods. Manually stock adjusting each medication (We wouldn't recommend this) as well as the legacy carry forward tool in the office section.

**Note:** Once toggled on the eMARx care home application will now automatically carry stock over to the new cycle. This will happen on the first sync of the new cycle. If the system fails to carry stock into the new cycle, then a support ticket is created. This feature can be switched off in the same way.



Go to Stock Control.



It will be switched off by default. Toggle on and sync the device.